

HIGH FLOW INSUFFLATION UNIT

Optimum Abdominal Conditions Allow Focused Surgery



Clear Interior Spaces For Clear-Thinking Procedures

The UHI-3 lets you keep your full concentration where it belongs. Its smoke evacuation maintains clear operative views, while the 35 L/min. insufflation quickly responds to gas leaks. In addition, the large front panel display allows you to check critical data at a glance.

Smoke Evacuation

Smoothly Removes Smoke & Mist


Just press the foot switch (optional) to evacuate smoke and mist generated during electrosurgery and ultrasonic procedures, so you can keep your eyes on your surgery.



Automatic Mist Evacuation Upon Ultrasonic Output

When the UHI-3 is connected to the SonoSurg Ultrasonic Surgical System, mist evacuation begins simultaneously with ultrasonic output. (Ultrasonic output is controlled by the MAJ- 51 SonoSurg-G Foot Switch for the SonoSurg Ultrasonic Surgical System.)





Warnings & Automatic Pressure Relief

Light And Sound Alarms Alert You To Problems. Excess Gas Is Automatically Released When Abdominal Pressure Becomes Too High.

Automatic Relief

When abdominal pressure becomes too high, a warning lamp and alarm are activated, and a pinch valve

on the suction tube opens automatically to release excess gas.



AWARNING . EXCESSIVE PRESSURE . TUBE OBSTRUCTED

Warning Alarms

Upon detecting a tube obstruction or closed stopcock, a warning lamp and sound alarm are activated. When the CO₂ gas cylinder is empty or the high pressure hose is improperly connected, the bar graph display for the gas supply shows a red light and an alarm sounds.

AWARNING . EXCESSIVE PRESSURE . TUBE OBSTRUCTED

Powerful 35 L/min. Insufflation


Responds To Abdominal Gas Leaks By Quickly Returning To The Preset Level.



The UHI-3 features a powerful 35 L/min. and high gas-supply. It can quickly respond to gas leaks and return to the preset rate.

• 3 Main Flow Rate Settings

You can quickly set flow rate to LOW 0.1-1.0, MED 1.1-19, or HIGH 20-35 L/min. 0.1-step settings between 0.1 and 2.0 L/min. allow delicate control of gas supply.



• Connectable To Medical Gas Pipeline

The UHI-3's unique pressure control mechanism lets you connect it to your hospital's medical gas pipeline, eliminating the need for gas cylinders.



Low-Noise Design

Noise is kept especially low through an innovative opening and closing mechanism in the electromagnetic valve.



Large Front Panel Display


Check Critical Data At A Glance

Actual and preset pressures are displayed simultaneously. You can easily monitor variations in abdominal pressure and other necessary information, while staying focused on what you're doing.





System Chart



Standard Set

- UHI-3 High Flow Insufflation Unit
- MAJ-591 Suction Tube

Includes Tube with Green Line, Small Diameter Tube (transparent), Differential Diameter Connector, Y-Shaped Connector, Extension Tube, and Small Luer-Lock Connector

- MAJ-590 Insufflation Tube Includes Tube (transparent) and Filter Connecting Tube
- Power Cord

Options



• MH-317 Foot Switch

- MAJ-1080 Cylinder Hose (PIN) for UHI-3
- Includes Hose for Cylinder, Cap, and Wrench MAJ-1081 Cylinder Hose (DIN) for UHI-3 Includes Hose for Cylinder, Cap1, Cap 2, Wrench (large) and Wrench (small)
- MAJ-1082 Cylinder Hose (ISO) for UHI-3
- Includes Hose for Cylinder, Cap1, Cap 2, Wrench (large), and Wrench (small) • MAJ-1084 Medical Gas Pipeline Adaptor (NIST) for UHI-3
- v/ Wrench MAJ-1085 Medical Gas Pipeline Adaptor (DISS) for UHI-3 w/ Wrench

Specifications

Item		Specifications
Applicable gas		CO2 gas for medical use (Connectable to CO2 cylinder or medical gas pipeline)
Abdominal pressure	Indicators	Digital displays; Bar graph
control	Settings	Abdominal pressure can be adjusted in 1 mmHg increments between 3 and 25 mmHg. To set pressure for 21 mmHg at 20 mmHg, press abdominal pressure control switch continuously for 3 seconds. Previous pressure setting is retained. Selected pressure setting is placed into memory. For safety reasons, settings at 21 mmHg or higher are memorized as 20 mmHg.
Flow rate settings	3 flow setting modes	*HIGH* mode: 20 to 35 L/min. *MED* mode: 1.1 to 19 L/min. *LOW* mode: 0.1 to 1.0 L/min. When power is turned ON, *LOW* mode is selected.
VOLUME indicator		Digital displays (0.0 to 19.9 L: 0.1 L steps / 20 to 999 L: 1 L steps) *Pressing the reset button will return the display to 0.0 L.
Power requirements	Voltage (AC)	100 - 240 V
	Frequency	50/60 Hz
	Input	150 VA
	Voltage fluctuation	Within ±10%
Classification (Electro-Medical Equipment)	Type of protection against electric shock	Class I (3-pin power cord)
	Degree of protection against electrical shock of applied part	Type CF Applied Part Insufflation tube (MAJ-590, MH-397) Suction tube (MAJ-591)
Size	Dimensions	295 (W) x 340 (D) x 150 (H) mm (housing dimensions)
	Weight	10 kg
Smoke evacuation function (optional)		When the foot switch (optional) is depressed, insufflation and suction are interlocked to maintain abdominal pressure at the set pressure while smoke evacuation is performed. Inoperative when "LOW" is selected as insufflation mode. Inoperative when the abdominal pressure is below 3 mmHg.

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.



OLYMPUS MEDICAL SYSTEMS CORP. OLYMPUS WINTER & IBE GMBH OLYMPUS MEDICAL SYSTEMS EUROPA GMBH OLYMPUS Store Failed Weiter Stars 14-16, 2009 Handlig, Gennary One Corporate Drive, Orangeburg, New York, Nover, Joseph Corporate Drive, Orangeburg, New York, Nover, Joseph Corporate Drive, Orangeburg, New York, New Yor

6100 Blue Lagoon Drive, Suite Stor Milann, Lister T-KEYMED LTD. Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

OLYMPUS SINGAPORE PTE LTD. OLYMPIC AND SHOLE SALES CONTROL TO THE SALES AND SALES A Dep 100 NOT Supervision And Instance And Ins OLYMPUS MOSCOW LIMITED LIABILITY COMPANY OLYMPUS AUSTRALIA PTY LTD

Printed in Japan R312UB-1005