



Tabletop steam sterilizer
for medical and dental professionals





NEW INTERFACE. EASY TO USE. ITALIAN DESIGNED.

A complete new experience thanks to the blue, backlit capacitive touchpad with acoustic signals to keep the user informed. The matt finish covers and scratch-free aluminium details give your practice a new look! Available in B & S Class and 18 or 23 litre chamber capacity, it's the perfect solution for high throughput of medical & dental instruments.

4 KG OF STERILE MATERIAL IN JUST 20 MINS

Hera offers 10 pre-set cycles to meet all your needs. Choosing the cycle according to the type of instruments to be sterilized allows you to save time, money and protect your instruments.

The water loading pump allows distilled water to be pumped from a tank placed on the floor directly to the autoclave.

All incoming water is tested to make sure it meets EN 13060 standard helping to prolong the life of your machine.

The Vacuum safety locking technology allows you to safely close the door without unnecessary efforts.

Opening it is even easier, just press 'Door Release'.

RELIABLE

Thanks to European components, the Newmed Hera autoclave lasts longer.

SAFE

The self-control system alerts you when maintenance is required.

INTELLIGENT

The sensors control the operating parameters and adapt the cycle time to internal and external conditions.

ADAPTABLE

Up to 10 preset, clinically validated sterilization programs.

FUNCTIONAL

Simple use and easy maintenance.

PERFECT FOR YOUR NEEDS

The chamber is made of thick AISI 304 steel without any welding to prevent leakage.

ACCURATE

The standard water quality sensor constantly monitors water quality. Good water quality guarantees perfect sterilization and longer life to your instruments.



FAST CYCLES. SLEEK DESIGN. RELIABLE PROCESSING.

Our sterilizers are innovative and incorporate leading Quiksteam™ technologies designed for high performance. The DentiSure & Kronos sterilizer range also have a simple and easy-to-use interface. Available in B, S or N Class and sterilising chamber capacity ranging from 6 to 60 litres, we are sure our autoclaves will meet your requirements.

SAVE TIME

Time is money for any business. Our intelligent software and patented QuikSteam™ Technology, enables you to process your instruments faster than ever before. With some of the fastest B class cycles on the market. The Autoclave has a useful start delay function, you can now process your instruments quickly and effectively.

SAVE MONEY

Your sterilizer is the heart of your practice, if it stops working, you stop working which effects your profitability. Midmark & Newmed autoclaves are renowned for outstanding reliability. Coupled with highly trained engineers who provide technical support, and a responsive customer support team who are there to support your every need, we will not let you down. Overhead costs can also affect your

business' bottom line. On our latest models, we improved efficiency by reducing water and electricity consumption by more than 20%, saving you money.

WE PUT YOU AT THE HEART OF OUR BUSINESS

Our sterilizers are the result of continuous feedback from our customers to our experienced research and development team, who are forever striving to provide you with the instrument processing products and services you deserve. Your trust is very important to us, therefore our units come as standard with all the accessories you would expect including the tray holder, trays, removal tool and water connections included. We ensure low cost of ownership by offering you a remarkably affordable quality maintenance programme



DentiSure B6



DentiSure B18 and B23



DentiSure B60SI

PLEASE VISIT OUR WEBSITE FOR PRODUCT INFORMATION

 newmed Hera



www.dentalexmedical.co.uk/hera

 midmark[™] DentiSure



www.dentalexmedical.co.uk/dentisure

 newmed Kronos



www.dentalexmedical.co.uk/kronos

CONTACT FORM



www.dentalexmedical.co.uk/contact-us

 midmark[™]

 newmed

UK Responsible Person MHRA

DENTALEX Ltd.

Singleton Court Business Centre

Wonastow Road West


Monmouth

NP25 5JA

Phone: 01600 719 919

Email: info@dentalexmedical.co.uk

Web: www.dentalexmedical.co.uk

 www.linkedin.com/company/dentalex-in-partnership-with-newmed

