

CE<sub>1023</sub>

Wireless Portable Dental X-Ray System

# BPD-I

*BEST QUALITY & BEST DESIGN – The best veterinary partner*



- Slim and light design for dental treatment*
- Wide and high resolution LCD*
- Easy Graphic Operation and Ergonomic Design*
- Low patient dosage and excellent image quality*
- Large capacity battery*
- Wide exposure ranges (0.01~2.00 sec.)*
- Compatible with chemical film and digital sensor*



## Wireless Portable Dental X-Ray System

# BPD-I



### Specifications

X-Ray Generator	DC High Frequency	
X-Ray Installation	Tube voltage	60 kV
	Tube current	2 mA
X-Ray Tube	Anode type	Stationary
	Anode angle	20°
	Focal spot size	0.8 mm
	Heat capacity	8.5 KHU
Exposure time	0.01~2.00 sec.	
Total Filtration	1.8 mm Al	
SSD	Min.100 mm	
Total Weight	2.3 Kg	
Dimension	174(W) X 143(H) X 150(D) mm	
Power requirement	Rechargeable Battery	
	22.2 VDC	