



# Digital Portable X-ray System

## Minibloc EPX-Series

- High Performance, High Frequency **EPX-Series** Portable X-ray system
- Compact & Lightweight
- Provides high resolution imaging at lower dosages
- Durable design with solid body
- Soft Touch Controls with Digital displays
- Calibrated cassette size indicator dials
- Available with Inverted control panel for table use
- Two stage dynamic, Auto line compensation
- Compatible with DR systems
- Dual Integrated Laser Pointer



## Digital Portable X-ray System



EPX-F1600



EPX-F2400



EPX-F2800



EPX-3200



EPX-F4000



EPX-F5000

### Specification

	EPX-F1600	EPX-F2400	EPX-F2800	EPX-3200	EPX-F4000	EPX-F5000
Max Output Power	1.6kW	2.4kW	2.8kW	3.2kW	4kW	5kW
Generator Frequency	100kHz	70kHz	70kHz	70kHz	70kHz	70kHz
kV Range	40~90kV	40~100kV	40~120kV	40~100kV	50~110kV	50~110kV
mA Range	12~30mA	16~40mA	12~40mA	25~60mA	20~80mA	20~100mA
mAs Range	0.4~100mAs	0.4~100mAs	0.4~100mAs	0.4~100mAs	0.4~100mAs	0.4~100mAs
Focal Spot	1.2mm*1.2mm	1.2mm*1.2mm	1.2mm*1.2mm	2.0mm*2.0mm	1.8mm*1.8mm	1.8mm*1.8mm
Target Angle	19°	16°	16°	16°	15°	15°
Anode Heat Storage	10kHU	50kHU	50kHU	40kHU	42kHU	42kHU
Input Voltage	110V/220V	110V/220V	110V/220V	200V~240V	200V~240V	200V~240V
Phase/Frequency	Single/50~60Hz	Single/50~60Hz	Single/50~60Hz	Single/50~60Hz	Single/50~60Hz	Single/50~60Hz
Dimension(mm)	335*153*170	360*200*190	360*200*190	360*200*190	360*200*190	360*200*190
Packing Size(mm)	520*400*330	520*400*330	520*400*330	520*400*330	520*400*330	520*400*330
Net weight	6.8kG	11kG	12kG	12kG	12kG	12kG
G-weight	20kG	22kG	22kG	22kG	22kG	22kG

### Mobile Stand For Portable X-ray Unit



EMS-100



EMS-100



EMS-H



EMS-H



Note: Ecotron reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.



Aquarius Centre, Edison Close, Waterwells Business Park,  
 Quedgeley, Gloucestershire, GL2 2FN  
 Tel: 01452 347101  
 Email: info@photonsurgicalsystems.co.uk  
 Web: www.photontonsurgicalsystems.co.uk