


# BeneHeart R3 Vet

Veterinary Electrocardiograph

## Your Faithful Consultant for Resting ECG Diagnosis




[www.mindrayanimal.com](http://www.mindrayanimal.com)  
[service@mindrayanimal.com](mailto:service@mindrayanimal.com)

## Reliable Results

- A 1,000Hz high digital sampling rate and 16,000Hz for pacemaker signal detection.
- Together with different filter setting, BeneHeart R3 Vet provides a reliable ECG waveform record comes with diagnostic results.


## Unique Recorder



Compatible with both rolling paper and Z-fold paper, you can easily switch between these two styles of papers without dismantling the pressure lever.


## Great Mobility

BeneHeart R3 Vet weighs only 1.2kg with battery, clinicians can easily carry the machine with single hand and perform resting ECG test in different scenarios.



## Long Operation Time

Large capacity rechargeable Lithium ion battery is provided, enough for 500 records printing or a 6-hour continues work without printing.



## Clear Display

5-inch color screen come with a high resolution of 800\*480 pixels, offering a clear view of ECG waveforms with detail recorded.

## Informatic Solution

Patient data can be saved as different formats, and transferred to PC or network server for historical data backup, review, printing and tele-diagnosis.

