型号	Veta 3S
Gas Supply	02
Flow Meter	Mechanical, 0-4LPM
Vaporizer Position	One position, ISO/SEV opt.
Vaporizer	V60 Vaporizer with Flow/Pressure /Temperature compensation function
Breathing System	Integrated
APL Valve	0-70 cmH2O, With one-button block
APL Valve ACGO	0-70 cmH2O, With one-button block Standard, One-button swtich
	·
ACGO	Standard, One-button swtich














www.mindrayanimal.com service@mindrayanimal.com

P/N:ENG-Veta 3S-210285X4P-20221209





Veta 3S Veterinary Anesthesia Machine

Safety and simplicity



Integrated Breathing System

- Anesthetic delivery faster and more stable with smaller system volume
- No tubes exposed, reducing the risk of leakage
- One step to remove Soda lime canister without any tools





High Quality Vaporizer

- Mindray V60 vaporizer, with flow /pressure /temperature compensation function
- Accurate output and good air tightness
 - Maintenance-free

APL Valve

- APL valve with quantitative scale
- One-button block, rising 30 cm H₂O





Easy to operate

 ACGO one-button switch between Non-rebreathing and rebreathing system

Anesthesia Gas Filter Canister

• Efficient adsorption of anesthesia waste gas to protect the safety of veterinarians and create a safe operating room





Excellent Movability

- A Four-wheeled trolley flexible and stable, easy to move
 - Caster with brake, stable and reliable